- Femur neck osteolysis.
- Insertion in growing cartilage.
- Secondary osteolysis of the segment.

Varus osteotomy plate

- Varus osteotomy.

9. Warnings

- Repeat use of a single use device is strictly forbidden. Each device used once shall be disposed of according to the hospital procedures. This is the same even where it appears to be intact; the device may have small faults or internal stresses that if the item is used may lead to fatigue failure.
- It is strictly prohibited to carry out any modification whatsoever on devices. If this recommendation is not followed, Citieffe disclaims any liability for any subsequent consequences.
- Any devices that have been taken out of its sterile package but have not been used or have become soiled should be handled according to hospital protocol. Citieffe does not recommend the reprocessing/resterilization of sterile or soiled devices.
- Do not use damaged, scratched, improperly treated devices or devices handled by unauthorized persons.
- Mixing of devices from different suppliers is not recommended for reasons of metallurgy, mechanics and design. Citieffe disclaims any liability in case of implants from different sources being mixed.

10. Suggestions

- Possible allergies or other reactions to the implant materials, even rare, should be considered, analysed (if necessary) and excluded before surgery.
- A preliminary review of the special surgical instruments is recommended.
- Use only titanium screws with titanium plates, and stainless steel screws with stainless steel plates.
- Particular attention should be paid to the right insertion of the blade when using varus osteotomy plate.
- Use only original Citieffe devices and accessories.
- The patient should be informed about the possibility of a new fracture even when the device is still implanted or just after its removal. This possibility decreases when the space, left empty, is filled again.
- At the end of the treatment and upon surgeon decision, once the devices have been removed, they should be disposed of as medical special waste.

11. Possible adverse effects

- Device bending or breakage due to an excessive intra or post-operative stress.
- Allergic reaction or rejection.
- Mobilization and loss of bone fixation.
- Discomfort, pain or paresthesia generated by implants.
- Nerve damage due to surgical trauma.
- Infection.
- Secondary fractures.

There are risks related to the failure of the implants. The surgeon and the patient must be aware of these risks and possible need of an early revision surgery.

12. MRI Safety Information

- The devices have not been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment. The devices have not been tested for heating, migration, or image artifact in the MR environment.
- The safety of the patient in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
- The patients should systematically report that they are wearer of implants. The patients should be informed of this recommendation.

13. Traceability

To ensure the traceability of the devices, the user should attach the labels provided in the packaging to the patient's clinical records.

14. Surgical technique

For the operative technique see the surgical technique manual of each product.

The Surgeon should make sure that the Patient reports immediately any adverse reactions, even if different from the ones indicated above.

For further information, please refer to Citieffe.

15. Symbols

	CATALOGUE NUMBER		USE BY (YEAR / MONTH / DAY)
	LOT NUMBER		SINGLE USE, DO NOT RE-USE
	ETHYLENE OXYDE STERILIZED		CAUTION: CONSULT INSTRUCTION FOR USE
	GAMMA RAYS STERILIZED		CONSULT INSTRUCTIONS FOR USE
	MANUFACTURER		DO NOT USE IF PACKAGE IS DAMAGED
	DATE OF MANUFACTURE (YEAR / MONTH / DAY)		NON STERILE

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