Citieffe srl a socio unico
Via Armadori, 21
40012 Calderara di Reno (Bologna) - Italy
Tel. +39 051 721850 - Fax +39 051 721870
info@citieffe.com - www.citieffe.com

ITALIANO

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INSTRUCTIONS FOR USE OF ESTREMO

Please, read before use.

This package insert does not include all of the information necessary for selection and use of a

device. Please see full labelling for all necessary information. Federal law (U.S.A.) restricts this device to be sale by or on the order of a licensed physician.

1. Description and intended use

The Nails Estremo - Citieffe Nailing System and Estremo - Ankle Fusion are single-use medical devices to be used as intramedullary internal fixation for the stabilization of long bone and tibiotalocalcaneal arthrodesis.

The Estremo Cortical Screws, Estremo Condylar Screws and Estremo Cancellous Screws are single-use medical devices to be used to lock the nail in the medullary canal and/or stabilize bone fragments.

The Estremo Locking End caps and Estremo Protection End caps are single-use medical devices to be used to protect the screwed extremity of the intramedullary nails and/or lock the screws. The range includes medical devices from Class Ib compliant with Regulation EU 2017/745. For further information regarding models, dimensions and sizes see the Citieffe catalogues.

2. Materials

The material used for devices manufacturing is Ti6Al4V ELI titanium alloy (ASTM F-136, ISO 5832-3 standards). The device material is provided on the packaging label.

3. Sterilization

The devices are provided sterile and have been exposed to validated gamma rays irradiation (R) sterilization process.

Sterile devices shall not be sterilized or resterilized by the user before use.

The manufacturer and the distributor accept no liability for sterilization of the devices performed by the user.

4. General safety and performance notes

- For the operative technique see the surgical technique manual of each product.
- The devices are designed according to state-of-the-art knowledge of implant technology and the applicable rules of engineering skills. However, the safety and functionality of devices can only be guaranteed if surgeons observe and comply with the notes in this instruction for use and the surgical technique.
- The sole purpose of devices is to support the healing process. They are no replacement for intact tissue and bone material. Due to the anatomy of human bones, implants are subject to limits with regard to their size and thickness.
- The tolerance, the resistance and the duration of the system depend on the material, on the correct use and assembling of the components and on the use of the dedicated instrument set.
- For further information and to request the summary of safety and clinical performance of the medical device, please refer to Citieffe.

5. Storage, handling and treatment of devices

Always check the expiration date (shelf life) and the indicator of sterility on the sterile packaging before opening. If the expiration date has past or if the indicator of sterility indicates non-sterile, the devices should be considered non-sterile and unusable.

- The sterile packaging should be checked for defects in the sterile barrier before opening. If any defects are noticed, the devices should be considered non-sterile and unusable.
- Open the sterile packaged device only during surgery, in the sterile field and immediately before use. Follow the hospital protocol in order to maintain the sterility of the device during handling. Remove all protective packaging before implantation.
- Visually inspect the devices before use. Any devices with damage, scratches, flaws, corrosion, discoloration, residue or debris should be discarded.
- Use extreme care in the handling of devices. Cutting, bending or scratching the surface of metal components can significantly reduce the strength and fatigue resistance of the device.
- Do not machine or modify devices in any way, unless required by the surgical technique and allowed by the manufacturer. In case of doubt, request a written recommendation from the manufacturer.
- Store products packed in their original packaging.
- Users should always wear appropriate personal protective equipment when handling the devices.
- Storage areas for devices should be clean from dust, away from areas of humidity, away from sunlight and at room temperature. This recommendation is equally valid for the transport and packaging of the devices.

6. Intended users

- The devices shall only be handled and used by qualified and trained healthcare practitioner (e.g. surgeons who are experienced in orthopedic surgery, operating room staff, and individuals involved in preparation of the device).
- Surgery shall only be performed by operating surgeons who are experienced in orthopedic surgery and who have acquired the appropriate qualifications, are aware of general risks of orthopedic surgery and are familiar with the product-specific surgical procedures.
- The surgeon's education, training and professional judgment must be relied upon to choose the most appropriate implant type and size and treatment option. No mandatory training on the product is required before using the Estremo implants.
- Citieffe cannot be held liable for complications arising from incorrect indication, surgical techniques or asepsis; this is the sole responsibility of the surgeon.

7. Target group of patients

The intended target patient group for internal synthesis treatment with Estremo consists of all skeletally mature patients who have a clinical condition attributable to one of the indicators of use for which the medical devices are intended.

From predinal evaluation, state of the art, and post market surveillance data, there are no known limitations of use related to the special elderly population (>65).

However, it is not recommended to use the device on patients who have one or more of the contraindications listed in these instructions for use.

Moreover, further precautions, especially in vulnerable populations (e.g. pregnant and breast-feeding populations) or suspected contraindicated population, should be taken as indicated in section "Suggestions" 11, keeping in mind also the risks outlined in section 12.

8. Indications

The Nail Estremo - Citieffe Nailing System is indicated for femoral fracture fixation, which may include the following: diaphyseal femoral fractures, supra-condylar femoral fractures, corrective osteotomies, pseudoarthrosis, non-unions and mal-unions.

The Nail Estremo - Citieffe Nailing System is indicated for tibia fracture fixation, which may include the following: diaphyseal tibia fractures, corrective osteotomies, pseudoarthrosis, non-unions and mal-unions.

The Nails Estremo - Citieffe Nailing System and Estremo - Ankle Fusion are indicated for tibiotalocalcaneal arthrodesis.

The Estremo Screws and End caps are indicated for use with Estremo Nails for fracture fixation and tibiotalocalcaneal arthrodesis.

9. Contraindications

- In progress infection or infection risk due to present or previous pathologies in the area to be treated.
- Hypersensitivity or reactivity or allergy to the material.
- Insertion through the growing cartilages.
- Estremo - Ankle Fusion: Femur fractures and tibia fractures.
- Conditions which tend to pre-empt the patient's ability or willingness to restrict activities or follow direction during treatment.
- Bone stock compromised by disease, infection or prior implantation that cannot provide adequate support and/or fixation of the devices.
- Obesity. An overweight or obese patient can produce loads on the implant which can lead to failure of the fixation of the device or to failure of the device itself.
- Any mental or neuromuscular disorder which could create an unacceptable risk of fixation failure or complications in postoperative care.
- Any drug addiction and/or tendency to abuse drugs, medicines or alcohol.
- Other medical or surgical conditions which would preclude the potential benefit of surgery.

10. Warnings

- The devices are not approved or intended for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.
- Repeat use of a single use device is strictly forbidden. Each device used once shall be disposed of according to the hospital procedures. This is the same even where it appears

to be intact; the device may have small faults or internal stresses that if the item is used may lead to fatigue failure.

- It is strictly prohibited to carry out any modification whatsoever on devices. If this recommendation is not followed, Citieffe disclaims any liability for any subsequent consequences.
- It is strictly forbidden to recondition or resterilize sterile devices taken out of the original packaging. Any devices that have been taken out of its sterile package but have not been used or have become soiled should be handled according to hospital protocol.
- Do not use damaged, scratched, improperly treated devices or devices handled by unauthorized persons.
- Mixing of devices from different suppliers is not recommended for reasons of metallurgy, mechanics and design unless otherwise specified (see the operative technique). Citieffe disclaims any liability in case of implants from different sources being mixed.

11. Suggestions

- Possible allergies or other reactions to the implant materials, even rare, should be considered, analysed (if necessary) and excluded before surgery.
- A preliminary review of the special surgical instruments is recommended.
- Use only titanium screws and end caps with titanium nails.
- Use only original Citieffe devices and accessories.
- The devices are designed for temporary implantation until the bone consolidation occurs. Therefore, if no bone consolidation occurs or if consolidation is not sufficient, the system may fail.
- The Devices are not intended for full weight bearing in patients with complex unstable fractures until sufficient bone consolidation is confirmed in the follow up X-rays.
- The patient should be informed about the possibility of a new fracture even when the device is still implanted or just after its removal. This possibility decreases when the space, left empty, is filled again.
- At the end of the treatment and upon surgeon decision, once the devices have been removed, they should be disposed of as medical special waste. In case of no protusion or pain, it is not recommended to remove the implants in elderly and weaken people.

12. Residual risks, adverse events and possible undesirable side effects

- Device bending or breakage due to an excessive intra or post-operative stress.
 - Incarceration or out coming of screws and protection end caps.
 - Allergic reaction or rejection.
 - Discomfort, pain or parestesia generated by implants.
 - Nerve damage due to surgical trauma.
 - Infection.
 - Limb length inequality caused by implants poor positioning or choice.
 - Bone damage, including secondary fractures, osteolysis or bone necrosis.
- There are intrinsic risks associated with anesthesia and surgery and residual risks related to the failure of the implants or complications. The latter risks may be caused by improper use of the device with respect to this manual and the intended surgical technique, use in the presence of contraindications, incorrect device selection and placement, reuse of a disposable device. The surgeon and the patient must be aware of these risks and possible need of an early revision surgery.

Any incidents associated with the use of the product, including those other than those mentioned above, shall be reported without undue delay to the manufacturer and the national authority.

Manufacturer's contact for reporting:

- support@citieffe.com

13. Patient information

- The patient must be informed of the residual risks, adverse events and possible undesirable side effects associated with the use of the medical device, as well as the implications for carrying capacity, mobility and general living conditions.
- The patient must be informed of the possibility of a refracture occurring even when the implant is in place or immediately after its removal; the possibility is less later, when the space left free by the implant is filled.

14. MRI Safety Information

INFORMATION NOT VALID FOR US MARKET: Non-clinical testing has demonstrated that titanium intramedullary nails can be scanned safely after implantation of the device under the following conditions:

- Static magnetic field: 3T;
- Maximum spatial field gradient: 153 T/m/;
- Maximum value of the specific absorption rate (SAR) reported by the MR system as a whole-body average: 2W/Kg.

Under the scan conditions defined above, titanium intramedullary nails produce a temperature rise of 0,3 °C at 3 T after 15 minutes of continuous scanning.

Non-clinical testing has demonstrated that titanium intramedullary nails scanned under the conditions defined above with static magnetic field of 3T do not under significant displacement forces and torque effects and produce image artifacts up to 38.26 mm from the implant.

For further information regarding the safety and compatibility of the devkes in the MR (Magnetic Resonance) environment, the Surgeon should contact Citieffe.

15. Traceability

To ensure the traceability of the devices, the user should attach the labels provided in the packaging to the patient's clinical records.

16. Symbols

Symbols in the label are harmonized and in accordance with EN ISO 15223-1.

Non-harmonized symbols:

R_X Only "Caution: Federal law (USA) restricts this device to sale by or on the order of a physician"

The symbols glossary is available at this address <https://www.citieffe.com/symbolsglossary>

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