

Modular monolateral limb reconstruction system




Modular and versatile:
various configurations possible

The compression/distraction mechanism
is an integral part of the fixator body:
simple control during treatment

Stainless steel, self-tapping bone screws
with double diameter thread
for a long-lasting stability

Small line (fixator body and clamps)
indicated for pediatric use

*Simple solutions
of complex reconstructive problems*



Indications

Rekrea external fixators are indicated in adult or pediatric patients in reconstructive procedures for the treatment of:

- shortening of the upper and lower limbs;
- deformity correction or bone-loss with the techniques of lengthening;
- bone transport and compression/distraction at one or multiple levels.

- ✔ Integrated compression/distraction unit, no additional items
- ✔ Independent movement of individual clamps
- ✔ All components packaged sterile and single-use
- ✔ For adult and pediatric patients



Clinical case

Biphocal bone transport

Courtesy of Dr. Carlos Satizabal
Clínica de Trauma y Reconstrucción
Ósea, Hospital Militar Central
(Bogotá)



Application of the fixator
with 4 clamps



Double osteotomy



Biphocal bone
transport