



Sansone

Revision Acetabular Cup



prosthesis



pelvis

citieffe[®]
Essential moves in Trauma

Sansone | Revision Acetabular Cup



prosthesis



pelvis

Easy to position and anchor in the iliac bone

Orientable HA coated
polar screw with angular stability



Optimal bone anchoring independent
of the cup position



BR-017A2.EN - Rev. 0 - 02/2016



One or more peripheral screws



Rotational control



Safety ring



Guaranteed locking of the polar screw

Indications

In case of bone defects due to prosthetic mobilization and severe congenital and / or acquired secondary morphological alterations, for example:

- outcomes of pelvic trauma
- mobilization of the acetabulum with bone defect (according to the AIR classification: type II, III, IV)
- congenital hip dysplasia
- severe rheumatic arthropathy

citieffe[®]
Essential moves in Trauma



Citieffe s.r.l.

Via Armaroli, 21 - 40012 Calderara di Reno (BO) - Italy
Tel +39 051 721850 - Fax +39 051 721870
info@citieffe.com - www.citieffe.com

CE
0476