

ORDERING INFORMATION

PK-130 Elbow Fixator Procedure Kit

STERILE 

CONTENTS



1 x **F4-2260**

Elbow Articulated External Fixator

Contents:

1 articulated external fixator (pin holders $\phi 6$ mm)

1 K. wire $\phi 2 \times 200$ mm

1 locking pin

NOTE: unique fixator for right and left side



1 x **F4-2749**

Elbow Fixator R.O.M. limiter



4 x **F4-140120** Cortical bone screw $\phi 4$ mm L120-25 mm

4 x **F4-150160** Cortical bone screw $\phi 5$ mm L160-35 mm


Code	Thread ϕ mm	Stem ϕ mm	Thread L. mm	Total L. mm
F4-140120	4	5	25	120
F4-150160	5	5	35	160












1 x **SF1381**

Targeting guide for K. wire $\phi 2$ mm

ST.A.R.90 F4 Elbow Articulated External Fixator



	2 x SF1397	Trocar, short
	2 x SF1396	Trocar, long
	2 x SF1392	Cannula, short (80 mm)
	2 x SF1391	Cannula, long (110 mm)
	2 x SF1334	Drill bit ø3 mm (AO joint)
	2 x SF1335	Drill bit ø4 mm (AO joint)
	1 x SF1305	Quick-connect adapter for bone screw (stem ø5mm)
	1 x SF1356	Handle hexagon 6 mm
	1 x DF900005	Handle ø5 mm